

**Product Selection Guide**

<u>Products</u>	<u>Absorbency</u>	<u>Size</u>	<u>Dimensions / Inches (cm)</u>	<u>ML</u>	<u>OZ</u>	<u>Case Qty</u>	<u>Reorder #</u>
<b>MoliMed Classic</b> 		Mini	Length x Width ~11x5 Inches (27x13.5cm)	140 ml	4.7 oz	18/14	168634
		Midi	Length x Width ~13x6 Inches (34x16cm)	212 ml	7.2 oz	12/14	168644
		Maxi	Length x Width ~18x8 Inches (45x20cm)	404 ml	13.7 oz	12/14	168654
<b>MoliForm Premium</b> 		Normal	Length x Width (OAL) ~24x13 Inches (62x32.5cm)	550 ml	18.6 oz	4/30	168019
		Plus	Length x Width (OAL) ~24x13 Inches (62x32.5cm)	720 ml	24.3 oz	4/30	168219
		Extra	Length x Width (OAL) ~27x14 Inches (69x36cm)	960 ml	32.5 oz	4/30	168319
		Super	Length x Width (OAL) ~27x14 Inches (69x36cm)	1360 ml	46.0 oz	4/30	168919
<b>MoliPants 2000</b> 	N/A	Medium	Waistband Color Blue = Medium	N/A	N/A	8/25	947791
	N/A	Large	Waistband Color Brown = Large	N/A	N/A	8/25	947792
<b>MoliCare Extra</b> 		X-Small X-S	Waist Range ~16-28 Inches (40-70cm)	860 ml	29.0 oz	4/30	169248
		Small Size 1	Waist Range ~20-32 Inches (50-80cm)	1040 ml	35.2 oz	3/30	169448
		Medium Size 2	Waist Range ~28-47 Inches (70-120cm)	1240 ml	41.9 oz	3/30	169648
		Large Size 3	Waist Range ~40-59 Inches (100-150cm)	1430 ml	48.3 oz	3/30	169848
		X-Large Size 4	Waist Range ~55-69 Inches (140-175cm)	1600 ml	54.0 oz	4/14	169948
<b>MoliCare Super</b> 		Small Size 1	Waist Range 20-32" (50-80cm)	1360 ml	46.0 oz	3/30	169450
		Medium Size 2	Waist Range 28-47" (70-120cm)	1700 ml	57.5 oz	3/30	169650
		Large Size 3	Waist Range 39-59" (100-150cm)	1870 ml	63.2 oz	3/30	169850

**Product Features Guide**

	<i>Soft Cloth Backsheet</i>	<i>Nonwoven Topsheet</i>	<i>Anti-Leak Guards</i>	<i>Leg Elastics</i>	<i>Blue Barrier</i>	<i>Wetness Indicator</i>	<i>Sizing Lines</i>	<i>Fastening Tapes, 2-Piece</i>	<i>Absorbent Core, 3-Part</i>	<i>Adhesive Strip</i>	<i>Dry-Plus Layer</i>	<i>High-Dry SAP with CyDex</i>	<i>Compact Packaging</i>	<i>Dermatological-Tested</i>	<i>Latex-Free</i>
<b>MoliMed Classic</b>	▲	▲		▲						▲	▲	▲	▲		▲
<b>MoliForm Premium</b>	▲	▲	▲	▲	▲	▲			▲		▲	▲	▲	▲	▲
<b>MoliPants 2000</b>	Used with MoliMed and MoliForm * For multiple use * Wash @ 60 degrees C or 140 degrees F * Colored waistband for side identification Long leg sections for added comfort * Soft, fine and strong elastic material for secure fit and comfort														
<b>MoliCare Extra</b>	▲	▲	▲	▲		▲	▲	▲	▲		▲	▲	▲	▲	▲
<b>MoliCare Super</b>	▲	▲	▲	▲		▲	▲	▲	▲		▲	▲	▲	▲	▲



**Quality Life Services, Inc.**

980 Nairn Ave  
Winnipeg, Manitoba R2L 0Y2  
Telephone: (204) 663-3562 Toll Free: (800) 363-0769  
Fax: (204) 663-1439 Toll Free Fax: (866) 613-6376  
Web Site: [www.qualitylifecares.com](http://www.qualitylifecares.com)

12-770 Industrial Avenue  
Ottawa, Ontario K1G 4H3  
Telephone: (613) 247-7554 Toll Free: (800) 731-6899  
Fax: (613) 247-9315 Toll Free Fax: (877) 229-2259  
Web Site: [www.qualitylifecares.com](http://www.qualitylifecares.com)